

Wheaton Sterile 1 Liter PET Media Bottle

NEW PRODUCT *Splash*

New Improved Design Features Wider Mouth, Heavier Construction, and a Double Seal Closure System

Thicker Wall Weight

The container wall weight has been increased by 30%, providing added handling strength and increased permeation barrier.

Wider Mouth Design

The mouth of the container has been increased from 33mm to 38mm, providing greater space for pipette entry, minimizing the risk of contamination.

Leak-Resistant Linerless Closure System

We integrated our superior leak-proof designed, double-seal closure system, currently on our HDPE leak-resistant bottles, into the 1000mL PET Media Square. The improved closure system greatly increases seal integrity and eliminates the risk of contamination by cap lining materials.

For bottle dimensions and quantity per case, see catalog #W219980 below.

Sterile PET Square Media Bottles

(Including new 1L)

- Permanent in-mold graduations
- No-drip pour lip
- Lightweight
- 20% headspace for additives

Wheaton PET Square Media Bottles offer the clarity of glass, the strength of plastic and a space-saving design. These sterile bottles are tested superior for pH stability, temperature durability, cloning efficiency and cytotoxicity. The bottles come with caps pre-attached and are shrink-wrapped in convenient trays.

Ordering Information

| Cat. No. | Size (mL) | L x W x H (mm) | Screw Cap Size | Qty/Tray | Qty/Case |
|----------|-----------|----------------|----------------|----------|----------|
| 219975 | 125 | 53 x 53 x 99 | 33-430 | 24 | 48 |
| 219979 | 500 | 75 x 75 x 175 | 33-430 | 12 | 24 |
| W219980 | 1000 | 91 x 91 x 220 | 38-430 | 12 | 24 |



WHEATON

Science Products

a division of WHEATON Industries Inc.

1501 North 10th Street, Millville, NJ 08332-2038

1.800.225.1437 ■ 1.856.825.1100

www.wheatonsci.com

Sterile 1 Liter PET Media Bottles are trademarks of Wheaton Industries Inc.
©Copyright, Wheaton Industries Inc., 2007 Printed in USA

Lit.No. 7796 07/07