

Socorex Profiller™ *electro* 446 Pipette Controller

NEW PRODUCT

Splash

- Powerful pump – fills 25 mL in 3 seconds
- Simple, fast operation
- Comfortable finger triggers
- Low battery light
- NiMH battery rechargeable during use
- Interchangeable protection membrane filter
- Stand suitable for bench or wall mounting
- One year warranty

This pipette controller fits any hand in a very comfortable way. The shape is ideal for working with glass or plastic pipettes within laminar flow hoods or microbiological safety cabinet. A powerful, yet quiet motor with a switch to set aspiration and dispensing speeds enables maximal working efficiency even with large pipetting volumes.

The variable speed switch enables selection of fast or slow aspiration and dispensing (with blow-out). Fast speed filling up to 25 mL in under 3 seconds. The battery charger enables charging battery either in storage position or during use. The convenient stand is suitable for bench or wall mounting with or without a pipette in place.

Ordering Information

Supplied with charger, two spare hydrophobic filters (0.45 and 0.2 µm), bench/wall stand and instruction manual.

Cat. No.	Description	Qty/Case
W851438	Profiller™ <i>electro</i> 446	1
W851436	Spare bench stand for Profiller™ <i>electro</i> 446	1
WI050950	Repl. 1.2V NiMH battery	1
W851444	Fast charger and cord	1
Replacements		
W835022	Membrane filter, 0.45µm	5
W851454	Membrane filter, 0.2µm	5
W835024	Silicone pipette holder	1
W835026	Nose cone, white	3
W835028	Alternative nose cones pink, violet, green	1 ea



WHEATON

Science Products

a division of WHEATON Industries Inc.

1501 North 10th Street, Millville, NJ 08332-2038

1.800.225.1437 ■ 1.856.825.1100

www.wheatonsci.com

Profiller is a trademark of Socorex ISBA
©Copyright, Wheaton Industries Inc., 2007 Printed in USA

