

NEW PRODUCT

Splash

Wheaton PET Diagnostic Bottles

- Crystal clear lightweight sterilized plastic containers
- Break resistant alternative to glass
- Ideal for reagents and buffer solutions
- Good resistance to some solvents, weak acids and oil
- Clear and amber available with caps attached
- Sterile polypropylene caps offered separately

Wheaton PET diagnostic bottles and caps are gamma radiation sterilized. These bottles are ideal for storage and shipment of reagents and buffer solutions. They are lightweight and break resistant. They come with pre-attached white polypropylene caps lined with F217 foamed LDPE. The bottles are packaged 20 per module and 5 modules per case. Sterile, colored caps can also be purchased separately. Assorted caps come with 1 pack of 20 of each color, blue, red, yellow, green, and black in a case of 100.

| Cat. No. | Size (mL) | Color | Qty/Pk | Qty/Cs |
|--------------|-----------|-------|--------|--------|
| Clear | | | | |
| W219981 | 5 | Clear | 20 | 100 |
| W219982 | 10 | Clear | 20 | 100 |
| W219983 | 20 | Clear | 20 | 100 |
| Amber | | | | |
| W219985 | 5 | Amber | 20 | 100 |
| W219986 | 10 | Amber | 20 | 100 |
| W219987 | 20 | Amber | 20 | 100 |

Polypropylene Screw Caps

| Cat. No. | Cap Size (mm) | Color | Qty/Pk | Qty/Cs |
|----------|---------------|----------|--------|--------|
| W219990 | 20-410 | Blue | 20 | 100 |
| W219991 | 20-410 | Red | 20 | 100 |
| W219992 | 20-410 | Yellow | 20 | 100 |
| W219993 | 20-410 | Green | 20 | 100 |
| W219994 | 20-410 | Black | 20 | 100 |
| W219995 | 20-410 | White | 20 | 100 |
| W219996 | 20-410 | Assorted | 20 | 100 |



WHEATON

Science International

a division of WHEATON Industries Inc.

1501 North 10th Street, Millville, NJ 08332-2038 USA

+1.856.825.1100 ■ +1.856.825.1368 (F)

www.wheatonsci.com