

**Coming Soon...**

**NEW PRODUCT**

*Splash*

# Wheaton RediMix Cap™

*Excellent for medical, diagnostic and lyophilization applications*

- Dispenses and mixes controlled amount of diluent without removing the cap
- Ensures complete and precise mixing without a pipette
- Saves users time and money

The new and revolutionary RediMix Cap™ provides the ability to store a controlled amount of diluent in the cap until you are ready to mix it with the product in the vial. The cap can hold between 1 – 2.0mL of diluent and can dispense and mix the diluent without having to remove the cap from the vial. This reduces the number of exposures to contaminants and ensures accuracy during mixing.

RediMix Cap™ is the answer which makes reconstitution simple, effective and convenient. It is great for use in diagnostic kit applications that require mixing as a step in the preparation of the final product. The ability to delay mixing increases shelf life.

Pipetting a measured amount of diluent leaves opportunity for human error. The RediMix Cap™ removes the need for a pipette and a second vial of diluent saving time and money while eliminating waste. Order the RediMix Cap™ pre-loaded with Water for Injection, custom diluents, or empty so you can load it with your own diluent.

**Capacity:** 1-2.0mL

**Container Style:** Fits any size vial or bottle with a 20 mm serum finish

For more information contact your local Wheaton Science Products Sales Representative.



← Vial not included



**WHEATON**

Science Products

a division of WHEATON Industries Inc.

1501 North 10th Street, Millville, NJ 08332-2038

1.800.225.1437 ■ 1.856.825.1100

[www.wheatonsci.com](http://www.wheatonsci.com)