

Buy 3 Cases of WHEATON Sample Vials and Receive 1 FREE Case*

Sample Vials in Lab File™

- Vials with caps attached packaged in partitioned trays
- Clear and amber vials made from borosilicate glass
- Screw caps designed for use with aqueous or organic samples

Vials for Aqueous Samples (with caps attached)

Cat. No.	Size (mL)	Size (dr)	Dia x H (mm) (cap on)	Cap Size	Qty/Case
Clear Vials with Solid Top Black Phenolic 14B Rubber Lined Caps					
224881	2	1/2	12 x 38	8-425	288
224882	4	1	15 x 48	13-425	144
224884	8	2	17 x 63	15-425	144
224885	12	3	19 x 68	15-425	144
224886	16	4	21 x 73	18-400	144
224888	24	6	23 x 88	20-400	144
225288	20	5	28 x 60	24-400	72
225289	25	6.25	28 x 73	24-400	72
225290	40	10	28 x 98	24-400	72

Amber Vials with Solid Top Black Phenolic 14B Rubber Lined Caps

224981	2	1/2	12 x 38	8-425	288
224982	4	1	15 x 48	13-425	144
224984	8	2	17 x 63	15-425	144
W224820	20	5	28 x 60	24-400	72
W224840	40	10	28 x 98	24-400	72

Vials for Organic Samples (with caps attached)

Cat. No.	Size (mL)	Size (dr)	Dia x H (mm) (cap on)	Cap Size	Qty/Case
Clear Vials with Solid Top Black Phenolic Teflon® Faced 14B Rubber Lined Caps					
W224581	2	1/2	12 x 38	8-425	288
W224582	4	1	15 x 48	13-425	144
W224584	8	2	17 x 63	15-425	144
W224585	12	3	19 x 68	15-425	144
W224586	16	4	21 x 73	18-400	144
W224588	24	6	23 x 88	20-400	144
W224589	20	5	28 x 60	24-400	72
W224590	25	6 1/4	28 x 73	24-400	72
W224591	40	10	28 x 98	24-400	72

Amber Vials with Solid Top Black Phenolic Teflon® Faced 14B Rubber Lined Caps

W224681	2	1/2	12 x 38	8-425	288
W224682	4	1	15 x 48	13-425	144
W224684	8	2	17 x 63	15-425	144
W224604	20	5	28 x 57	24-400	72
W224605	40	10	28 x 95	24-400	72

LIMITED TIME OFFER

Redemption Coupon

Please send my **FREE Case of Sample Vials Catalog No.** _____.
Please enclose or fax a copy of your proof of purchase: (invoice, packing slip, purchase order).

Name: _____

Company: _____

Street: _____

City: _____ State: _____


Zip/Postal: _____ Country: _____

Tel: _____

E-Mail: _____

This offer is void where prohibited by law or company policy. All information must be completed to be a valid entry. Offer valid in North America only thru 04/30/10. Offer limited to two redemptions per customer. *Free case must be of equal or lesser value.

Lit. No. 8390 01/10



Sample Vials in Lab File™

Offer Valid Thru April 30, 2010

To place an order, contact your local laboratory supply distributor.
To receive your **FREE Case of Sample Vials**, complete the redemption coupon and send or fax/email it along with your proof of purchase to:

Wheaton Science Products
Attention: Fulfillment Coordinator
1501 North 10th Street ■ Millville, NJ 08332-2038 USA
Email: jackie.garrison@wheaton.com
1.800.225.1437 (T) ■ 1.856.825.9035 (F)