



WHEATON

CRIMPenstein™

Electric Crimper/Decapper



Advanced Ergonomic Design
Engineered to Save Time While Making
Your Crimping and Decapping Easier



Lightweight and easy to operate

Enables handheld semi-automatic operation for extended periods with less effort.

Quick change jaw-sets available in 8 mm, 11 mm, 13 mm, and 20 mm sizes

Provides cost savings by eliminating the need to purchase separate crimpers and decappers.

Electric powered for quick setup and high volume operations

Eliminates the need for regulated compressed air or batteries which discharge over time.

Adjustable jaw intensity with memory

Allows user to regulate jaw action as required to achieve the perfect crimp.

Electrify your crimping and decapping operations with the new innovative Wheaton Crimpenstein. Engineered to improve productivity and efficiency in a variety of laboratories performing chromatography, QC testing, lyophilization, vaccine development or any other applications where vials are capped with aluminum seals. This AC powered crimper/decapper provides fast, consistent results with a simple push of a button. Jaw action intensity can be adjusted as needed. The ability to change jaw-sets offers an unmatched level of flexibility and economy for applications where various sizes of aluminum seals are required. A sturdy stand is available which will hold the Crimpenstein in a fixed vertical position for operations where vials need to be handled one at a time. The new Flip Cap Crimper Jaw attaches 13 and 20 mm flip cap aluminum seals. Flip cap seals are easily removed using the standard 13 or 20 mm Decapper Jaw.

Crimpenstein Crimper and Decapper



Crimpenstein with Crimper and Decapper Jaw-Set

- Includes 1 Crimpenstein Electric Crimper/Decapper Unit, 1 Crimper Jaw, 1 Decapper Jaw, 1 Power Cord, Quick Start Card and Instruction Manual
- Packaged in convenient storage/carrying case

| Part No. | For Use With | Size (mm) | Voltage | Qty/Case |
|-----------|-----------------------------|-----------|---------|----------|
| W225600-A | Standard Aluminum Seals | 8 | 120 VAC | 1 |
| W225601-A | Standard Aluminum Seals | 11 | 120 VAC | 1 |
| W225602-A | Standard Aluminum Seals | 13 | 120 VAC | 1 |
| W225603-A | Standard Aluminum Seals | 20 | 120 VAC | 1 |
| W225643-A | NEW Flip Cap Aluminum Seals | 13 | 120 VAC | 1 |
| W225644-A | NEW Flip Cap Aluminum Seals | 20 | 120 VAC | 1 |

Crimpenstein with Crimper Jaw

- Includes 1 Crimpenstein Electric Crimper/Decapper Unit, 1 Crimper Jaw, 1 Power Cord, Quick Start Card and Instruction Manual
- Packaged in convenient storage/carrying case

| Part No. | For Use With | Size (mm) | Voltage | Qty/Case |
|-----------|-----------------------------|-----------|---------|----------|
| W225605-A | Standard Aluminum Seals | 8 | 120 VAC | 1 |
| W225606-A | Standard Aluminum Seals | 11 | 120 VAC | 1 |
| W225607-A | Standard Aluminum Seals | 13 | 120 VAC | 1 |
| W225608-A | Standard Aluminum Seals | 20 | 120 VAC | 1 |
| W225637-A | NEW Flip Cap Aluminum Seals | 13 | 120 VAC | 1 |
| W225638-A | NEW Flip Cap Aluminum Seals | 20 | 120 VAC | 1 |



Storage/Carrying Case

Crimpenstein with Decapper Jaw

- Includes 1 Crimpenstein Electric Crimper/Decapper Unit, 1 Decapper Jaw, 1 Power Cord, Quick Start Card and Instruction Manual
- Packaged in convenient storage/carrying case

| Part No. | For Use With | Size (mm) | Voltage | Qty/Case |
|-----------|------------------------------------|-----------|---------|----------|
| W225610-A | Standard & Flip Cap Aluminum Seals | 8 | 120 VAC | 1 |
| W225611-A | Standard & Flip Cap Aluminum Seals | 11 | 120 VAC | 1 |
| W225612-A | Standard & Flip Cap Aluminum Seals | 13 | 120 VAC | 1 |
| W225613-A | Standard & Flip Cap Aluminum Seals | 20 | 120 VAC | 1 |



Crimpenstein Stand

- Holds the Crimpenstein Unit in vertical position for operations where vials need to be handled one at a time.
- Stores two Crimper or Decapper Jaws of any size in the base of stand

| Part No. | Description | Qty/Case |
|----------|--------------------|----------|
| W225630 | Crimpenstein Stand | 1 |

Crimpenstein Crimper and Decapper Jaws



Crimper and Decapper Jaw-Set, 20 mm
W225618

Crimper and Decapper Jaw-Set

- Expands capability of Crimpenstein Unit by adding different size Crimper and Decapper jaws that can be purchased separately
- Crimper Jaw attaches aluminum seal to bottle/vial and Decapper Jaw removes the aluminum seal
- Jaw-Set includes 1 Crimper and 1 Decapper Jaw of the same size

| Part No. | For Use With | Size (mm) | Qty/Case |
|----------|-------------------------|-----------|----------|
| W225615 | Standard Aluminum Seals | 8 | 1 |
| W225616 | Standard Aluminum Seals | 11 | 1 |
| W225617 | Standard Aluminum Seals | 13 | 1 |
| W225618 | Standard Aluminum Seals | 20 | 1 |
| W225641 | Flip Cap Aluminum Seals | 13 | 1 |
| W225642 | Flip Cap Aluminum Seals | 20 | 1 |



Crimper Jaw, 13 mm
W225622

Flip Cap Crimper Jaw, 20 mm
W225640

Crimper Jaw

- Expands capability of Crimpenstein Unit by adding different size Crimper Jaws
- Crimper Jaw attaches aluminum seal to bottle/vial
- Includes 1 Crimper Jaw

| Part No. | For Use With | Size (mm) | Qty/Case |
|----------|-------------------------|-----------|----------|
| W225620 | Standard Aluminum Seals | 8 | 1 |
| W225621 | Standard Aluminum Seals | 11 | 1 |
| W225622 | Standard Aluminum Seals | 13 | 1 |
| W225623 | Standard Aluminum Seals | 20 | 1 |
| W225639 | Flip Cap Aluminum Seals | 13 | 1 |
| W225640 | Flip Cap Aluminum Seals | 20 | 1 |



Decapper Jaw, 20 mm
W225628

Decapper Jaw

- Expands capability of Crimpenstein Unit by adding different size Decapper Jaw
- Decapper Jaw removes the aluminum seal from the bottle/vial, and can be used for both Standard and Flip Cap Aluminum Seals
- Includes 1 Decapper Jaw

| Part No. | For Use With | Size (mm) | Qty/Case |
|----------|------------------------------------|-----------|----------|
| W225625 | Standard & Flip Cap Aluminum Seals | 8 | 1 |
| W225626 | Standard & Flip Cap Aluminum Seals | 11 | 1 |
| W225627 | Standard & Flip Cap Aluminum Seals | 13 | 1 |
| W225628 | Standard & Flip Cap Aluminum Seals | 20 | 1 |

For a complete list of Serum Bottles, Vials and Closures, see our on-line catalog at www.wheatonsci.com



Crimpenstein Electric Crimper / Decapper - International Ordering Information

| Country Plug Code | North America (A) | Japan (B) | Continental Europe (C) | United Kingdom (D) | Australia/China (F) | Italy/Chile (G) | India (J) | Qty/Case |
|-------------------|-------------------|-----------|------------------------|--------------------|---------------------|-----------------|-----------|----------|
| Description | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | |

Crimpenstein for Standard Aluminum Seals

| | | | | | | | | |
|--|-----------|-----------|-----------|-----------|-----------|-----------|-----------|---|
| With 8 mm Crimper and Decrimper Jaw-Set | W225600-A | W225600-B | W225600-C | W225600-D | W225600-F | W225600-G | W225600-J | 1 |
| With 11 mm Crimper and Decrimper Jaw-Set | W225601-A | W225601-B | W225601-C | W225601-D | W225601-F | W225601-G | W225601-J | 1 |
| With 13 mm Crimper and Decrimper Jaw-Set | W225602-A | W225602-B | W225602-C | W225602-D | W225602-F | W225602-G | W225602-J | 1 |
| With 20 mm Crimper and Decrimper Jaw-Set | W225603-A | W225603-B | W225603-C | W225603-D | W225603-F | W225603-G | W225603-J | 1 |

Crimpenstein for Flip Cap Aluminum Seals

| | | | | | | | | |
|--|-----------|-----------|-----------|-----------|-----------|-----------|-----------|---|
| With 13 mm Crimper and Decrimper Jaw-Set | W225643-A | W225643-B | W225643-C | W225643-D | W225643-F | W225643-G | W225643-J | 1 |
| With 20 mm Crimper and Decrimper Jaw-Set | W225644-A | W225644-B | W225644-C | W225644-D | W225644-F | W225644-G | W225644-J | 1 |

Crimpenstein for Standard Aluminum Seals

| | | | | | | | | |
|------------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|---|
| With 8 mm Crimper Jaw | W225605-A | W225605-B | W225605-C | W225605-D | W225605-F | W225605-G | W225605-J | 1 |
| With 11 mm Crimper Jaw | W225606-A | W225606-B | W225606-C | W225606-D | W225606-F | W225606-G | W225606-J | 1 |
| With 13 mm Crimper Jaw | W225607-A | W225607-B | W225607-C | W225607-D | W225607-F | W225607-G | W225607-J | 1 |
| With 20 mm Crimper Jaw | W225608-A | W225608-B | W225608-C | W225608-D | W225608-F | W225608-G | W225608-J | 1 |

Crimpenstein for Flip Cap Aluminum Seals

| | | | | | | | | |
|------------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|---|
| With 13 mm Crimper Jaw | W225637-A | W225637-B | W225637-C | W225637-D | W225637-F | W225637-G | W225637-J | 1 |
| With 20 mm Crimper Jaw | W225638-A | W225638-B | W225638-C | W225638-D | W225638-F | W225638-G | W225638-J | 1 |

Crimpenstein for Standard & Flip Cap Aluminum Seals

| | | | | | | | | |
|--------------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|---|
| With 8 mm Decrimper Jaw | W225610-A | W225610-B | W225610-C | W225610-D | W225610-F | W225610-G | W225610-J | 1 |
| With 11 mm Decrimper Jaw | W225611-A | W225611-B | W225611-C | W225611-D | W225611-F | W225611-G | W225611-J | 1 |
| With 13 mm Decrimper Jaw | W225612-A | W225612-B | W225612-C | W225612-D | W225612-F | W225612-G | W225612-J | 1 |
| With 20 mm Decrimper Jaw | W225613-A | W225613-B | W225613-C | W225613-D | W225613-F | W225613-G | W225613-J | 1 |

For a complete list of Chromatography Vials and Closures, visit www.wheatonsci.com

