

WHEATON

Purchase 2 Cases of Crimp Top ValuePaks and Receive a 11mm E-Z Crimper™ FREE*



W225164



W225165



W225166



W225167



ValuePak Autosampler Vials, Crimp Top

- Vials are packaged 100 per shelf pack and 10 packs per case
- Manufactured from certified clear and amber Type I borosilicate glass
- 40% larger opening improves sample accessibility while reducing autosampler needle damage

Cat. No.	Size (mL)	Description	Qty/Case
E-Z Vial® with Aluminum Crimp Top Seal with PTFE/Red Rubber Liner			
W225164	1.8	Clear	1000
W225165	1.8	Clear with Writing Patch	1000
W225166	1.8	Amber	1000
W225167	1.8	Amber with Writing Patch	1000

E-Z Crimper™

- For attaching aluminum seals to bottles and vials with a crimp/serum finish
- Cushioned ergonomic handle reduces hand fatigue
- Labeled for quick E-Z Crimper size identification
- Polished crimping jaws provide consistent sealing
- 2 year warranty

Cat. No.	Description	Qty/Case
W225301*	Attaches 11mm Standard Aluminum Seals	1

*The 11mm E-Z Crimper is supplied with a plunger insert to give the option of crimping the top of a vial creating an extra tight seal. We recommend this only for standard opening vials.

Complete the redemption coupon below and send or mail/fax/email it along with your proof of purchase to:

Limited Time Offer Offer Valid thru June 30, 2012

WHEATON

Attention: Joel Torres
1501 North 10th Street
Millville, NJ 08332-2038 USA

Email: joel.torres@wheaton.com
800.225.1437 (T) | 856.825.9035 (F)

Please send my FREE 11mm E-Z Crimper
Cat. No. _____.

Please enclose or fax a copy of your proof of purchase: Invoice, packing slip, or purchase order with your redemption.

Name: _____
Company: _____
Address: _____
Address: _____
City: _____
State: _____
Zip/Postal: _____
Country: _____
Tel: _____
E-Mail: _____

This offer is void where prohibited by law or company policy. All information must be completed to be a valid entry. Offer valid in North America only thru 06/30/2012. Limited to two redemptions per customer. *Free case must be of equal or lesser value. Not for resale.

www.wheaton.com