

WHEATON

Buy 2 Calibrex™ 520 Bottle Top Dispensers, Get 1 FREE*



The Calibrex™ digital 520

Designed for performance handling of liquids from a large variety of bottles and flasks, the Calibrex dispensers combine the latest in dosing technology, high-tech materials and ergonomic design. They are manufactured and tested to fully comply with current safety regulations. Three color-coded 520-models with adjustable volumes are ranging 0.25 - 2mL, 1 - 5mL and 1 - 10mL.

Product Advantages

- Superior chemical resistance
- Long lasting performance stability
- Simple construction - no tool for disassembling
- Instant volume setting
- In-lab calibration
- Autoclavable at 121°C/250°F fully assembled
- Two-year warranty

Performance and Ordering Information - Calibrex™ 520

Calibrex 520 dispensers have a 32mm base thread. They come with three adapters (28, 40, 45mm) and feed tube 300mm.

Cat. No.	Volume	Division
844000	0.25 - 2mL	0.05mL
844002	1 - 5mL	0.1mL
844004	1 - 10mL	0.25mL

Cat. No.	Inaccuracy (E%)			Imprecision (CV%)		
	Min. Vol.	Mid. Vol.	Max. Vol.	Min. Vol.	Mid. Vol.	Max. Vol.
844000	< ± 3.0 %	< ± 1.8 %	< ± 0.6 %	< 0.5 %	< 0.35 %	< 0.1 %
844002	< ± 2.0 %	< ± 1.3 %	< ± 0.6 %	< 0.5 %	< 0.35 %	< 0.1 %
844004	< ± 1.5 %	< ± 1.1 %	< ± 0.6 %	< 0.5 %	< 0.35 %	< 0.1 %

* Specify desired fixed volume within instrument range when ordering. Performance values obtained with bi-dist. water at constant temperature (± 0.5°C) comprised between 20 and 25°C in accordance with ISO 8655.

Complete the redemption coupon below and send or mail/fax/email it along with your proof of purchase to:

Limited Time Offer
Offer Valid thru June 30, 2012

WHEATON

Attention: Joel Torres
1501 North 10th Street
Millville, NJ 08332-2038 USA

Email: joel.torres@wheaton.com
800.225.1437 (T) | 856.825.9035 (F)

Please send my FREE Calibrex™ Dispenser,
Cat. No. _____.

Please enclose or fax a copy of your proof of purchase: Invoice, packing slip, or purchase order with your redemption.

Name: _____
Company: _____
Address: _____
Address: _____
City: _____
State: _____
Zip/Postal: _____
Country: _____
Tel: _____
E-Mail: _____

This offer is void where prohibited by law or company policy. All information must be completed to be a valid entry. Offer valid in North America only thru 06/30/2012. Limited to two redemptions per customer. *Not for resale.

www.wheaton.com