

WHEATON

Buy 2 Cases of BOD Bottles, and get a FREE BOD Bottle Rack*



BOD Bottles

- Ideal for incubating diluted samples of sewage, sewage effluents, polluted waters, and industrial wastes to determine the amount of oxygen required during the stabilization of the decomposable organic matter by aerobic biochemical action
- We recommend the 300mL capacity bottles for the five-day BOD test as referenced in Standard Methods for the Examination of Water and Wastewater, 21st Edition, 5210B, 2005
- Specially designed bottle shoulder radius that sweeps all air from inside the bottle during filling
- The interchangeable stoppers have a tapered bottom that prevents air entrapment
- The bottles have a flared mouth to form a water seal around the stopper that prevents air from being drawn into the bottle during incubation
- Stopper joint is compatible with the probes of the leading meters for BOD and dissolved oxygen
- Large, permanently screened-on writing patch on all bottles
- Manufactured from USP Type I borosilicate glass

Cat. No.	Number Sequence	Barcoding	Dia x H (mm)	Stopper	Qty/Case
60mL Capacity					
227494-00G	Un-numbered	No	43 x 115	Glass Pennyhead	36
227494-00	Un-numbered	No	43 x 115	Glass Robotic	36
227494-01G	01-36	Yes	43 x 115	Glass Pennyhead	36
227494-01	01-36	Yes	43 x 115	Glass Robotic	36
227494-02G	37-72	Yes	43 x 115	Glass Pennyhead	36
227494-02	37-72	Yes	43 x 115	Glass Robotic	36
227494-03G	73-108	Yes	43 x 115	Glass Pennyhead	36
300mL Capacity					
227498	Un-numbered	No	69 x 143	Without Stopper	24
227497-00G	Un-numbered	No	69 x 165	Glass Pennyhead	24
227497-00	Un-numbered	No	69 x 165	Glass Robotic	24
227497-01G	01-24	Yes	69 x 165	Glass Pennyhead	24
227497-01	01-24	Yes	69 x 165	Glass Robotic	24
227497-02G	25-48	Yes	69 x 165	Glass Pennyhead	24
227497-02	25-48	Yes	69 x 165	Glass Robotic	24
227497-03G	49-72	Yes	69 x 165	Glass Pennyhead	24
227497-03	49-72	Yes	69 x 165	Glass Robotic	24
227497-04G	73-96	Yes	69 x 165	Glass Pennyhead	24
227497-04	73-96	Yes	69 x 165	Glass Robotic	24
227497-05G	97-120	Yes	69 x 165	Glass Pennyhead	24
227497-05	97-120	Yes	69 x 165	Glass Robotic	24

* For additional number sequences or specials contact WHEATON Customer Service.

BOD Bottle Rack

- Conveniently and safely transports BOD bottles
- Racks are stackable

Cat. No.	For Bottle Size	Holds	L x W x H (in)*	Qty/Case
89174-478	60mL	20 Bottles	13 x 10 x 7 in	1
89174-476	300mL	12 Bottles	13 x 10 x 7 in	1

*33.03 x 25.4 x 17.78cm

Complete the redemption coupon below and send or mail/fax/email it along with your proof of purchase to:

Limited Time Offer
Offer Valid thru June 30, 2012

WHEATON

Attention: Joel Torres
1501 North 10th Street
Millville, NJ 08332-2038 USA

Email: joel.torres@wheaton.com
800.225.1437 (T) | 856.825.9035 (F)

Please send my FREE BOD Bottle Rack,
Cat. No. _____.

Please enclose or fax a copy of your proof of purchase: Invoice, packing slip, or purchase order with your redemption.

Name: _____
Company: _____
Address: _____
Address: _____
City: _____
State: _____
Zip/Postal: _____
Country: _____
Tel: _____
E-Mail: _____

This offer is void where prohibited by law or company policy. All information must be completed to be a valid entry. Offer valid in North America only thru 06/30/2012. Limited to two redemptions per customer. *Not for resale.

www.wheaton.com