



Flat Bottom Culture Tubes

Close-Out Sale!

Purchase
While
Supplies
Last!



Culture Tubes & Caps sold separately

Close-Out Sale on Culture Tubes Only

Part Number: W012710

Description: 16.5 x 125 mm Flat Bottom Clear Glass
Screw Thread Culture Tube with a *15-415 Finish

Packaging:

- Corrugated case of 288 pcs/ two trays per carton
- Two shrinkwrapped modules of 144 pieces, packed finish up
- Each module arranged 9 rows of 16 inside a chipboard partition

List Price: \$60.34/case

Close-Out Price: \$37.25/case

***Specifications:**

15-415 finish:

"T" Dim	0.566" - 0.581" (14.38 - 14.76 mm)
"E" Dim	0.506" - 0.521" (12.85 - 13.23 mm)
"S" Dim	0.022" - 0.047" (0.56 - 1.19 mm)
"H" Dim	0.535" - 0.560" (13.59 - 14.22 mm)
Overall Length:	125.0 mm +/- 1.0 mm
Tube OD:	16.0 - 16.50 mm
Wall Thickness:	1.10 mm average

Screw Caps

Culture Tubes sold without caps. If a cap is needed, you can purchase our screw cap W000035, which is not included in the Close-Out Sale.

Part Number: W000035

Description: Black Phenolic Screw Cap with 14B Rubber
Liner with a 15-415 Finish

Packaging: 7,400 caps/case

Price: \$670.00/case

To place an order contact
Wheaton Science Products
Customer Service:

800-225-1437
Credit card orders only



WHEATON

Science Products

a division of WHEATON Industries Inc.

1501 North 10th Street, Millville, NJ 08332-2038
1.800.225.1437 • 1.856.825.1100 • 1.856.825.1368 (F)
www.wheatonsci.com

No changes or substitutions allowed. Offer valid while supplies last. Wheaton has the right to discontinue this offer at any time without notification. All Sales are final.

SUPERIOR PRODUCTS

PROVEN QUALITY

EXCEPTIONAL SUPPORT